



# PRODUCT CATALOG



## FLOOR TYPE FREEZE DRYER BFFT-317-B

## FLOOR TYPE FREEZE DRYER BFFT-317-B

### FREEZE DRYER

Compact, stainless steel unit with casters for easy movement, Floor type freeze dryer provides unmatched process accuracy and reliability in performance. It optimizes your space requirements and reduces energy consumption.



## SPECIFICATIONS

Model	BFFT-317-B
Condenser Temperature	-110°C
Freeze Drying area	0.36 m <sup>2</sup>
Bulk capacity	3.6 L
Tray	φ240mm, 8pcs
Round bottom flask	1000ml, 500ml, 250ml, 100ml, 50ml
Alt Name	Manifold type Electric heating shelf

## APPLICATIONS

Perishable Substances, Waste Products, Laboratory, Research, Proteins, Organic Tissues, Waste Products, Plant Material, Polymers, Pharmaceuticals, Nutraceuticals, Plant material.



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)