



PURE WATER SUPPLY SYSTEM BFG-301

PURE WATER SUPPLY SYSTEM BFG-301

RO WATER SYSTEM WATER PURIFIER



SPECIFICATIONS

Model	BFG-301
Description	Reverse Osmosis System
Mounting	Horizontal or Wall-mounted
Filters	Sediment Filter, Precarbon Filter, R/O membrane Filter
Requirement to water source	Municipal tap water
Principle	Reverse Osmosis
Quality of the purified water	<20µs/cm
Ability	7-10 Litre/Hour
Product size	440mmx400mmx140mm
Package size	500mmx460mmx200mm
Voltage	220V 50Hz
Weight	10 Kgs
Alt Name	RO Water system Water Purifier



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com