



## WATER DISTILLER BFG-401

# WATER DISTILLER BFG-401

WATER DISTILLER



## SPECIFICATIONS

Model	BFG-401
Power	750W
Voltage	220V/50Hz
Capacity	1.5L/h
Chamber size	780mm*200mm
Package size	290mm*290mm*390mm
Product size	245mm*285mm*510mm
N.W.	3.5kgs
G.W.	5.0kgs
Alt Name	Water Distiller



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)