



PLATE SEALER BFG-501

PLATE SEALER BFG-501

SEALER



SPECIFICATIONS

Model	BFG-501
Voltage	200-240V, 100-120V
Frequency	50/60Hz
Power	500W
Max sealing length	110mm
Maximum sealing width	250mm
Net weight	5Kg
Dimensions	430mmx340mmx250mm
Sealing Temperature	180
Alt Name	Sealer



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com