



## MEDICAL AUTOCLAVE BFG-901

# MEDICAL AUTOCLAVE BFG-901

FAST DRY AUTOCLAVE



Double load:

1. More instruments can be loaded in the S/S baskets randomly.
2. To make your work easier, more efficient and costs saved.

Perfect dry:

1. Feature TSD technology, the drying results are perfect.
2. Don't worry about the WET packages any more.

Super fast:

1. The fastest cycle with the full load only need 39 minutes.
2. To speed up the instruments reprocessing.

## SPECIFICATIONS

Model	BFG-901
Capacity	23L
Input	220±22V / 50HZ
Power	2000 W
Chamber size	25.4*45 Cm
Product size	61*47*47 Cm(W*L*H)
Installation size	60Cm
N.W.	58Kgs
G.W.	61Kgs
Output	Data Printer, Label Printer (Optional)
Alt Name	Fast Dry Autoclave

## FEATURES

Double load:

1. More instruments can be loaded in the S/S baskets randomly.
2. To make your work easier, more efficient and costs saved.




Perfect dry:

1. Feature TSD technology, the drying results are perfect.
2. Don't worry about the WET packages any more.

Super fast:

1. The fastest cycle with the full load only need 39 minutes.
2. To speed up the instruments reprocessing.





**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada  
Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)