



DENTAL AUTOCLAVE BFG1B4 (BADT-104)

DENTAL AUTOCLAVE BFG1B4

AUTOCLAVE



- 1.Safety Valve
Prevent chamber from over pressure
- 2.Bacteria Filter
High quality,2µm precision
- 3.Steam Generator
Patented design, high efficiency and reliable
- 4.Air Switch
Shut down the machine in case of electricity leakage or short circuit
- 5.Double Fans
Bring faster cooling and more efficiency sterilization
- 6.Water Steam Separator
Patented design, prevent vacuum pump from oil damaging

SPECIFICATIONS

| | |
|-------------------|---------------------|
| Model | BFG1B4 |
| Old Model | BADT-104 |
| Voltage | 220V / 50Hz |
| Power | 1750W |
| Capacity | 22L |
| Chamber size | Φ25*45cm |
| Product size | 70*44.5*39cm |
| Installation size | 60cm |
| N.W. | 50kgs |
| Printer | optional |
| Alt Name | Table-top autoclave |



FEATURES



- 1.Safety Valve
Prevent chamber from over pressure



2.Bacteria Filter
High quality, 2µm precision



3.Steam Generator
Patented design, high efficiency and reliable



4.Air Switch
Shut down the machine in case of electricity leakage or short circuit



5.Double Fans
Bring faster cooling and more efficiency sterilization





6. Water Steam Separator

Patented design, prevent vacuum pump from oil damaging



7. Brass Piping

Withstand high temperature and pressure safe and reliable



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com