



**PRE AND POST VACUUM CLASS B AUTOCLAVE
BFG1C1 (BAPV-201)**

PRE AND POST VACUUM CLASS B AUTOCLAVE BFG1C1

PRE AND POST VACUUM CLASS B AUTOCLAVE



- 1.Safety Valve
Prevent chamber from over pressure
- 2.Bacteria Filter
High quality,2µm precision
- 3.Steam Generator
Patented design, high efficiency and reliable
- 4.Air Switch
Shut down the machine in case of electricity leakage or short circuit
- 5.Double Fans
Bring faster cooling and more efficiency sterilization
- 6.Water Steam Separator
Patented design, prevent vacuum pump from oil damaging

SPECIFICATIONS

Model	BFG1C1
Old Model	BAPV-201
Description / V/Hz	220/50
Voltage / W	2500
Capacity / L	8
Chamber size / cm	Φ 17x32
Product size / cm	53.5x39.5x32
Installation Size / cm	50
Weight / kg	45
Alt Name	Pre and Post Vacuum Class B Autoclave



FEATURES



- 1.Safety Valve
Prevent chamber from over pressure



- 2.Bacteria Filter
High quality,2µm precision



3. Steam Generator
Patented design, high efficiency and reliable



4. Air Switch
Shut down the machine in case of electricity leakage or short circuit



5. Double Fans
Bring faster cooling and more efficiency sterilization



6. Water Steam Separator
Patented design, prevent vacuum pump from oil damaging



7. Brass Piping
Withstand high temperature and pressure safe and reliable



8. German vacuum pump (optional)
More powerful and less noise.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com