



PRE AND POST VACUUM CLASS B AUTOCLAVE BFG1C3 (BAPV-203)

PRE AND POST VACUUM CLASS B AUTOCLAVE BFG1C3

PRE AND POST VACUUM CLASS B AUTOCLAVE



- 1.Safety Valve
Prevent chamber from over pressure
- 2.Bacteria Filter
High quality,2µm precision
- 3.Steam Generator
Patented design, high efficiency and reliable
- 4.Air Switch
Shut down the machine in case of electricity leakage or short circuit
- 5.Double Fans
Bring faster cooling and more efficiency sterilization
- 6.Water Steam Separator
Patented design, prevent vacuum pump from oil damaging

SPECIFICATIONS

Model	BFG1C3
Old Model	BAPV-203
Description / V/Hz	220/50
Voltage / W	2500
Capacity / L	18
Chamber size / cm	Φ 25x35
Product size / cm	61x50x43
Installation Size / cm	80
Weight / kg	53
Alt Name	Pre and Post Vacuum Class B Autoclave



FEATURES



- 1.Safety Valve
Prevent chamber from over pressure



- 2.Bacteria Filter
High quality,2µm precision



3.Steam Generator

Patented design, high efficiency and reliable



4.Air Switch

Shut down the machine in case of electricity leakage or short circuit



5.Double Fans

Bring faster cooling and more efficiency sterilization



6.Water Steam Separator

Patented design, prevent vacuum pump from oil damaging



7.Brass Piping

Withstand high temperature and pressure safe and reliable



8. German vacuum pump (optional)

More powerful and less noise.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com