



## PLATE SEALER BFG1F1

# PLATE SEALER BFG1F1

SEALER



## SPECIFICATIONS

|                       |                    |
|-----------------------|--------------------|
| Model                 | BFG1F1             |
| Voltage               | 200-240V, 100-120V |
| Frequency             | 50/60Hz            |
| Power                 | 500W               |
| Max sealing length    | 110mm              |
| Maximum sealing width | 250mm              |
| Net weight            | 5Kg                |
| Dimensions            | 430mmx340mmx250mm  |
| Sealing Temperature   | 180                |
| Alt Name              | Sealer             |



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)