



ULTRASONIC CLEANER BFG1G1

ULTRASONIC CLEANER BFG1G1

ULTRASONIC CLEANER



SPECIFICATIONS

Model	BFG1G1
Power	75W
Voltage	220V/50Hz
Tank size	240mm*140mm*150mm
Capacity	5L
Temp. range	20-80°C
Time range	1~60min
N.W.	4.0kgs
G.W.	5.0kgs
Alt Name	Ultrasonic Cleaner



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com