



PRE AND POST VACUUM CLASS B AUTOCLAVE BFG1I1

PRE AND POST VACUUM CLASS B AUTOCLAVE BFG111

BASIC LINE AUTOCLAVE



Durable & Reliable:

a. Near 20 years of manufacture experience.

b. Complete quality control system.


Effectiveness:

Stable control of temperature and pressure to guarantee the effectiveness during heavy and continuous sterilizations.

SPECIFICATIONS

Model	BFG111
Capacity	23L
Input	220V±22V / 50HZ
Power	2300 VA
Chamber size	25.4*45 Cm
Product size	61*43.5*42.5 Cm(W*L*H)
Installation size	60Cm
N.W.	38Kgs
G.W.	40Kgs
Alt Name	Basic Line Autoclave

FEATURES



Durable & Reliable:

a. Near 20 years of manufacture experience.

b. Complete quality control system.

Effectiveness:

Stable control of temperature and pressure to guarantee the effectiveness during heavy and continuous sterilizations.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com