



FLOOR STANDING LOW SPEED REFRIGERATED CENTRIFUGE BFH-10201

LOW SPEED REFRIGERATED CENTRIFUGE BFH-10201

LARGEST CAPACITY REFRIGERATED CENTRIFUGE



1. Widely used in the field of blood bank, pharmaceutical factory and laboratory.
2. Brushless frequency motor, in great torque, free maintenance, no powder pollution, quick in speed up and down.
3. 7" LCD touch screen display which indicates the program, speed, time, RCF, temperature.
4. Micro computer control, there are 100 kinds of program and 10 kinds of acceleration and deceleration for your choice.
5. CFC-free refrigerant, pre-cooling function, double cycle cooling, cold and hot alternating easily, free environment pollution and precise in temperature control.
6. Electric lid lock, compact design, super speed and imbalance protection.
7. The centrifuge body is made of high-quality steel, safe and reliable.

SPECIFICATIONS

Model	BFH-10201
Max. Speed	8000rpm
Max. RCF	14551xg
Max. Capacity	8 x 2000 ml (16 x 500 ml blood bag)
Time Range	1-99 h 59 min 59 sec
RPM/RCF Convert	Yes
Noise (dB)	≤ 65
Temperature	-20-40°C
Acc/Dec	10 Kinds
Speed Accuracy	±20 r/min
Temperature Accuracy	±1°C
Voltage (V/Hz)	AC 220V/380V, 50Hz/60Hz
Dimensions (L x W x H)	940 x 890 x 1000 mm
Net Weight (Kg)	600 kg
Certificates	CE, ISO & Calibration report are available
Alt Name	Largest Capacity Refrigerated Centrifuge

ITEMS INCLUDED

No	Name	Max Speed (rpm)	Max Volume (ml)	Max. RCF (xg)
1	Swing Rotor	4000	8 x 2000 ml	5716

ACCESSORIES FOR PURCHASE

No	Name	Max Speed (rpm)	Max Volume (ml)	Max. RCF (xg)
1	Swing Rotor	4000	6 x 2400 ml	5330
2	Angle Rotor	8000	6 x 1000 ml	14551
3	Swing Rotor	4200	2 x 6 x 1000 ml	5680

FEATURES

1. Widely used in the field of blood bank, pharmaceutical factory and laboratory.
2. Brushless frequency motor, in great torque, free maintenance, no powder pollution, quick in speed up and down.
3. 7" LCD touch screen display which indicates the program, speed, time, RCF, temperature.

4. Micro computer control, there are 100 kinds of program and 10 kinds off acceleration and deceleration for your choice.
5. CFC-free refrigerant, pre-cooling function, double cycle cooling,cold and hot alternating easily, free environment pollution and precise in temperature control
6. Electric lid lock, compact design, super speed and imbalance protection.
7. The centrifuge body is made of high-quality steel, safe and reliable.
- 8.Equipped with an automatic lid lock,the wind shield lid is close & open automatically with door,no manual lifting,with anti-pinch function as well (optional).
9. Special programs for blood components separation Largest Capacity Refrigerated Centrifuge



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com