



BENCHTOP HIGH SPEED CENTRIFUGE BFH-10601

HIGH SPEED CENTRIFUGE BFH-10601

TABLE TOP HIGH AND LOW SPEED UNIVERSAL CENTRIFUGE



1. Widely used for laboratory, chemistry and hospital.
2. 10 kinds of acceleration and deceleration
3. More than 10 rotors for your choice
4. High speed, low noise, compact design and small vibration
5. User-friendly LCD display provides speed, time and temperature parameters. 100 Programs can be stored.
6. Brushless frequency motor which in great torque, no powder pollution, free maintenance, quick in speed up and down
7. Imbalance sensor and over speed sensor, the centrifuge body is made of high-quality steel, very safe.

SPECIFICATIONS

Model	BFH-10601
Max. Speed	18500rpm
Max. RCF	23990xg
Max. Capacity	4 x 600 ml
Time Range	0-99 h 99 min 59 sec
RPM/RCF Convert	Yes
Noise (dB)	≤ 65
Acc/Dcc	10 Kinds
Speed Accuracy	±20 r/min
Voltage (V/Hz)	AC 220V/110V, 50Hz/60Hz
Size (W x D x H)	685 x 500 x 385 mm
Net Weight (Kg)	75 kg
Certificates	CE, ISO & Calibration report are available
Alt Name	Table Top High and Low Speed Universal Centrifuge

ACCESSORIES FOR PURCHASE

No	Name	Max Speed (rpm)	Max Volume (ml)	Max. RCF (xg)
1	High speed centrifuge	18500 rpm	4 x 600 ml	23990 xg
2	Swing rotor	4000	4 x 600 ml (round cup)	3530
3	Swing rotor	4000	4 x 600 ml (square cup)	3450
4	Swing rotor	4000	4 x 7 x 50 ml	3470
5	Swing Rotor	4000	4 x 18 x 10 ml	3150
6	Swing Rotor	4000	4 x 30 x 7.5 ml	2840
7	Fixed Rotor	14000	4 x 8PCR	12070
8	Fixed Rotor	13000	6 x 8PCR	16080
9	Fixed Rotor	14000	8 x 8PCR	13390
10	Fixed Rotor	13000	12 x 8PCR	17220
11	Fixed Rotor	14000	40 x 0.5 ml	19960
12	Fixed Rotor	18500	12 x 1.5 ml	23990
13	Fixed Rotor	15000	24 x 1.5 ml	20600
14	Fixed Rotor	13500	30 x 1.5 ml	19340
15	Fixed Rotor	11000	48 x 1.5 ml	12840
16	Fixed Rotor	15000	16 x 5 ml	19350
17	Fixed Rotor	14000	6 x 10 ml	16460

18	Fixed Rotor	13000	12 x 10 ml	17030
19	Fixed Rotor	11000	16 x 10 ml	13950
20	Fixed Rotor	11000	8 x 15 ml	12730
21	Fixed Rotor	9000	12 x 15 ml	9590
22	Fixed Rotor	5000	24 x 15 ml	3500
23	Fixed Rotor	5000	30 x 15 ml	3830
24	Fixed Rotor	12000	6 x 30 ml	14000
25	Fixed Rotor	12000	6 x 50 ml (conical)	16050
26	Fixed Rotor	12000	6 x 50 ml (round)	15960
27	Fixed Rotor	11000	4 x 85 ml	13560
28	Fixed Rotor	10000	4 x 100 ml	10310
29	Fixed Rotor	8000	6 x 100 ml	7280
30	Fixed Rotor	4000	12 x 100 ml	2970
31	Fixed Rotor	12000	24 capillaries	15800
32	Swing Rotor	15000	4 x 5 ml	19920
33	Vertical rotor	16000	16 x 5 ml	16450
34	Vertical rotor	14000	8 x 30 ml	19750
35	Microplate rotor	4000	2 x 4 x 96 well	3210



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com