



MICRO CENTRIFUGE BFH-11201

MICRO CENTRIFUGE BFH-11201

MICRO CENTRIFUGE



1. This centrifuge is widely used in laboratories and hospitals.
2. This model can fit in 2ml, 1.5ml, 0.5ml and 0.2ml tubes.
3. The handheld centrifuge is very suitable for microtubule filtration and fast centrifugation, the micro blood cell separation, microbial sample processing, experimental partition centrifugal of PCR.
4. Light weight, small size and beautiful design.
5. Easy control, super speed protection.

SPECIFICATIONS

Model	BFH-11201
Max. Speed	12000rpm
Max. RCF	8550xg
Max. Capacity	
Electric Source Input	AC85-264V / 50/60Hz
Rate of Work Input	30W
Motor	DC24V
Control	PI / Power direct drive
Fuse	PPTC / self recovery fuse / No replacement needed
Rotor Capacity	2/1.5/0.5/0.2 ml x 8 + 0.2 ml x 8 x 4 (PCR)
Maximum Speed [r/min.]	12,000 r/min ±5%
Size [mm]: Length x Width x Height	158 x 172 x 123 mm
Noise at Highest Speed	≤47 dB(A)
Alt Name	Micro Centrifuge

ACCESSORIES FOR PURCHASE

No	Name	Capacity (ml x tube number)	Maximum Speed (r/min.)	Max RCF
1	Angle Rotor	8x1.5/2ml	12000rpm	8550xg
2		0.2mlx8x4(PCR)		
3	Adapters	0.5ml		
4	Adapters	0.2ml		



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com