



# PRODUCT CATALOG



## LABORATORY LOW SPEED REFRIGERATED CENTRIFUGE BFH-3801

# LABORATORY LOW SPEED REFRIGERATED CENTRIFUGE

## BFH-3801

TABLETOP LARGE CAPACITY CENTRIFUGE



1. Brushless DC frequency motor with simpler construction, more reliable performance, longer life and quietly running.
2. Automatically electric lid lock, super speed, over temperature protection and imbalance protection. The centrifuge body is made of high quality steel, safe and reliable.
3. Rotor is connected to spindle by specialized taper sleeve, loading simple and quick, no direction, safe.
4. Microprocessor control, digital display in the speed, time, RCF in operation, speed raising and reducing quick, operate simply.
5. 10 kinds of program stored in the memory, 10 kinds of accelerating and decelerating speed for your choice.
6. 3 layers protection of steel structure and get the ideal centrifugation result

## SPECIFICATIONS

Model	BFH-3801
Max. Speed	5000rpm
Max. RCF	4730xg
Max. Capacity	4x800ml
Time Range	0~99min59sec
Noise	≤ 58dBA
Speed Accuracy	±20rpm
Size (W x D x Hmm)	685x500x405mm
Net Weight(Kg)	67KG
Power supply	AC 220V 50HZ 10A
Package	Wooden box
Alt Name	Tabletop Large Capacity Centrifuge

## ACCESSORIES FOR PURCHASE

No	Name	Max speed (rpm)	Max Volume (ml)	Max RCF (g)
1	Swing out rotor	4000	4x800ml	3380
2	Swing out rotor	4000	4x500ml	3310
3	swing out rotor	5000	4x500ml	4730
4	swing out rotor	5000	4x100ml	3020
5	swing out rotor	4000	8x100ml	3020
6	swing out rotor	4000	8x50ml	3200
7	swing out rotor	4000	16x15ml	3200
8	swing out rotor	4000	24x15ml	3200
9	swing out rotor	4000	88x5ml vacuum tube	2480
10	swing out rotor	4000	160x5ml vacuum tube	2480
11	swing out rotor	4000	32x10ml vacuum tube	2790
12	swing out rotor	4000	32x7ml vacuum tube adapter	2790
13	swing out rotor	4000	32x5ml vacuum tube adapter	2790
14	swing out rotor	4000	48x10ml vacuum tube	2880
15	swing out rotor	4000	48x7ml vacuum tube	2760
16	swing out rotor	4000	48x5ml vacuum tube	2480

17	swing out rotor	4000	72x10ml vacuum tube	3140
18	swing out rotor	4000	88x7ml vacuum tube	3140
19	swing out rotor	4000	160x7ml vacuum tube	3140
20	Micro plate rotor	4000	4x4x96 well	2940
21	Micro plate rotor	4000	2x4x96 well	3210

## FEATURES

1. Brushless DC frequency motor with simpler construction, more reliable performance, longer life and quietly running.
2. Automatically electric lid lock, super speed, over temperature protection and imbalance protection. The centrifuge body is made of high quality steel, safe and reliable.
3. Rotor is connected to spindle by specialized taper sleeve, loading simple and quick, no direction, safe.
4. Microprocessor control, digital display in the speed, time, RCF in operation, speed raising and reducing quick, operate simply.
5. 10 kinds of program stored in the memory, 10 kinds of accelerating and decelerating speed for your choice.
6. 3 layers protection of steel structure and get the ideal centrifugation result.
7. It is consistent with CE security certification and ISO9001 quality certification.



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)