



LABORATORY LOW SPEED CENTRIFUGE BFH1BD1 (BCBL-306)

LABORATORY LOW SPEED CENTRIFUGE BFH1BD1

UNIVERSAL CENTRIFUGE



1. Microprocessor control, LCD color display with touch key screen in speed, time, RCF in operation, speed raising and reducing quick, operate simply.
2. Brushless DC frequency motor with simpler construction, more reliable performance, longer life and quietly running.
3. Electric lid lock, super speed, over temperature protection and imbalance protection.
4. Adapters are available by experiment requirements.
5. The centrifuge body is made of high quality steel and the centrifuge is with three layers protection of steel structure, safe and reliable.
6. CE, ISO 9001 & ISO 13485 approved.

SPECIFICATIONS

Model	BFH1BD1
Old Model	BCBL-306
Max. Speed	5000rpm
Max. RCF	4730xg
Max. Capacity	4x250ml
Time Range	0~99min59sec
Noise	≤ 58dBA
Speed Accuracy	± 20r/min
Voltage(V/Hz)	AC 220V 50HZ 5A
Size (W x D x Hmm)	568x400x320mm
Net Weight(Kg)	35KG
Package	carton box
Alt Name	Universal centrifuge

ACCESSORIES FOR PURCHASE

No	Name	Max Speed (rpm)	Max Volume (ml)	Max. RCF (xg)
1	Swing out rotor	4000	4x250ml (round cup)	2880
2	Swing out rotor	4000	4x250ml (square cup)	2830
3	Swing out rotor	5000	4x50ml	4730
4	Swing out rotor	5000	4x100ml	4730
5	Swing out rotor	4000	8x50ml	2950
6	Swing out rotor	4000	16x15ml	2950
7	Swing out rotor	4000	24x15ml	2950
8	Swing out rotor	4000	48x5ml vacuum tube	2480
9	Swing out rotor	4000	48x7ml vacuum tube	2760
10	Swing out rotor	4000	40x10ml vacuum tube	2880

No	Rotor No.	Name	Max Speed (rpm)	Max Volume (ml)	Max. RCF (xg)
1	No.31377	Swing out rotor	4000	32x10ml vacuum tube	2580
2	No.31377	Swing out rotor	4000	32x7ml vacuum tube adapter	2580
3	No.31377	Swing out rotor	4000	32x5ml vacuum tube adapter	2580
4	No.31377	Swing out rotor	4000	24x10ml vacuum tube	2580

5	No.31377	Swing out rotor	4000	24x7ml vacuum tube adapter	2580
6	No.31377	Swing out rotor	4000	24x5ml vacuum tube adapter	2580
7	No.31377	Swing out rotor	4000	16x10ml vacuum tube	2580
8	No.31377	Swing out rotor	4000	16x7ml vacuum tube adapter	2580
9	No.31377	Swing out rotor	4000	16x5ml vacuum tube adapter	2580
10	No.31376	Micro plate rotor	4000	2X3x48 well	2300
11	No.31376	Micro plate rotor	4000	2X2x96 well	2490



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com