




LABORATORY LOW SPEED CENTRIFUGE BFH1BH1 (BCBL-308)

LABORATORY LOW SPEED CENTRIFUGE BFH1BH1

CENTRIFUGE FOR GEL CARDS



1. Brushless DC motor to make sure free maintenance no powder pollution, micro computer control and LCD display the time, speed and RCF to let operation easily.
2. There is specialized step-ladder centrifugation technology; different step-ladders could be designed by different reagents.
3. There are several kinds of regular centrifugation programs; other centrifugation programs also could be designed by customers' demand.
4. Electronic lock, over speed and imbalance protection, there is a self-check protection program, safe and reliable, simple operation, scientific and practicable.
5. With CE & ISO9001 & ISO13485 Approval.

SPECIFICATIONS

Model	BFH1BH1
Old Model	BCBL-308
Max speed	3000rpm
Max RCF	1190xg
Max volume	24 cards/12 cards
Noise	≤55dBA
Timer	0~99min59s
Power supply	AC 220 50HZ 2A
Dimension (12 cards)	460*320*230mm
Net. Weight (12 cards)	17Kgs
Dimension (24 cards)	500*360*250mm
Net. Weight (24 cards)	19Kgs
Speed accuracy	±20rpm
Package	carton box
Alt Name	Centrifuge for Gel Cards

APPLICATIONS

It is widely used in the test of Blood Group Identification,Cross-Matching of Blood,Irregular Antibody Screening etc.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com