



GEL INSTRUMENT BFH1BI1

GEL INSTRUMENT BFH1BI1

BLOOD CARDS INCUBATOR



1. Keep 24-Gel Cards in temperature of 37 degree.
2. 2 independent incubation zones with a capacity for 12 Gel cards
3. Digital display to set up time and keep the temperature
4. Audible programmable alarm to notify end of incubation period

SPECIFICATIONS

Model	BFH1BI1
Capacity	24 Gel cards
Temperature	37°C
Timer	0-99min
Power	AC220V 50Hz / 110V 60Hz
Dimensions	380x220x115 mm
Alt Name	Blood Cards Incubator

APPLICATIONS

Gel Cards Incubator is matched with our Gel Cards Centrifuge to use and holds 2 independent incubation zones with a capacity for 12 Gel cards,. It can be suitable for many types of gel cards from different manufacturers.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com