



## LABORATORY HIGH SPEED REFRIGERATED CENTRIFUGE BFH1BN1 (BCBH-307)

# LABORATORY HIGH SPEED REFRIGERATED CENTRIFUGE

## BFH1BN1

### HEMATOCRIT CENTRIFUGE



1. Brushless DC motor, no carbon dust pollution, free maintenance.
2. Microprocessor control, small LCD display.
3. Automatic rotor identification function.
4. Electronic lock, over speed and imbalance protection.
5. 10 kinds of program stored in the memory, 10 kinds of accelerating and decelerating speed for your choice.
6. It is with CE & ISO9001 & ISO13485 Certificates.

## SPECIFICATIONS

|                    |                       |
|--------------------|-----------------------|
| Model              | BFH1BN1               |
| Old Model          | BCBH-307              |
| Max. Speed         | 12000rpm              |
| Max. RCF           | 15800xg               |
| Max. Volume        | 24 capillary tubes    |
| Time Range         | 0~99min59sec          |
| Noise              | ≤ 55dBA               |
| Speed Accuracy     | ±20rpm                |
| Size (W x D x Hmm) | 335*270*190mm         |
| Net Weight(Kg)     | 13KG                  |
| Power supply       | AC 220V 50HZ 2A       |
| Package            | Carton box            |
| Alt Name           | Hematocrit centrifuge |



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)