



LABORATORY HIGH SPEED REFRIGERATED CENTRIFUGE BFH1BW1 (BCBLR-402)

LABORATORY HIGH SPEED REFRIGERATED CENTRIFUGE

BFH1BW1

TABLE TOP HIGH SPEED REFRIGERATED UNIVERSAL CENTRIFUGE



1. Refrigerated centrifuge. Widely used for laboratory, chemistry and hospital.
2. Pre-cooling design, 10 kinds of acceleration and deceleration
3. More than 10 rotors for your choice
4. High speed, low noise, compact design and small vibration
5. User-friendly LED display provides speed, time and temperature parameters.
6. Brushless frequency motor which in great torque, no powder pollution, free maintenance, quick in speed up and down
7. Imbalance sensor and over speed sensor, the centrifuge body is made of high-quality steel, very safe.

SPECIFICATIONS

Model	BFH1BW1
Old Model	BCBLR-402
Max. Speed	18500rpm
Max. RCF	23990xg
Max. Volume	4x180ml
Time Range	0~99min59sec
RPM/RCF Convert	Yes
Noise	≤ 55dBA
Acc/Dec	10 Kinds
Temperature	-20°C-40°C
Speed Accuracy	±20rpm
Temperature Accuracy	±1°C
Net Weight(Kg)	110KG
Power supply	AC 220V/110V 50HZ/60HZ
Certificates	CE,ISO & Calibration report are available
Alt Name	Table Top High Speed Refrigerated Universal Centrifuge

ACCESSORIES FOR PURCHASE

No	Name	Max Speed (rpm)	Max Volume (ml)	Max RCF (xg)
1	Swing Rotor	4000	4x180ml	2580
2	Microplate rotor	4000	2x2x96 hole	2020
3	Swing rotor	5000	4x1x50ml	4420
4	Swing rotor	5000	4x1x100ml	4420
5	Swing rotor	4000	4x4x15ml 4x2x50ml	2950
6	Swing rotor	4000	4x6x10/7ml vacuum tube	2580
7	Swing rotor	4000	4x6x5ml vacuum tube adapter	2170
8	Swing rotor	4000	4x4x10/7ml vacuum tube	2580
9	Swing rotor	4000	4x4x5ml vacuum tube adapter	2170
10	Fixed Rotor	14000	4x8PCR	12070
11	Fixed Rotor	13000	6x8PCR	16080
12	Fixed Rotor	14000	8x8PCR	13390
13	Fixed Rotor	13000	12x8PCR	17220
14	Fixed Rotor	14000	40x0.5ml	19960

15	Fixed Rotor	18500	12x1.5ml	23390
16	Fixed Rotor	15000	24x1.5ml	20600
17	Fixed Rotor	13500	30x1.5ml	19340
18	Fixed Rotor	11000	48x1.5ml	12840
19	Fixed Rotor	15000	16x5ml	19350
20	Fixed Rotor	14000	6x10ml	16460
21	Fixed Rotor	13000	12x10ml	17030
22	Fixed Rotor	11000	16x10ml	13950
23	Fixed Rotor	11000	8x15ml	12730
24	Fixed Rotor	9000	12x15ml	9590
25	Fixed Rotor	5000	24x15ml	9590
26	Fixed Rotor	5000	30x15ml	13800
27	Fixed Rotor	12000	6x30ml	14400
28	Fixed Rotor	12000	6x50ml(falcon)	16050
29	Fixed Rotor	12000	6x50ml(round)	15960
30	Fixed Rotor	11000	4x85ml	13560
31	Fixed Rotor	10000	4x100ml	10310
32	Fixed Rotor	8000	6x100ml	7280
33	Swing Rotor	15000	4x5ml	19920
34	Vertical rotor	16000	16x5ml	16450
35	Vertical rotor	14000	8x30ml	19750



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com