



MICRO CENTRIFUGE BFH1BY1 (BCBH-304)

MICRO CENTRIFUGE BFH1BY1

MICRO HIGH SPEED CENTRIFUGE



1. Brushless motor, no pollution, free-maintenance.
2. Microprocessor control, LCD display which indicates the speed, time, RCF in operation, 10 kinds of brake setting, operate simply.
3. Electric lid lock, super speed, imbalance protection. The centrifuge body is made of high quality steel, stainless steel chamber, safe and reliable.
4. Rotor is connected to spindle by specialized taper sleeve, loading simple and quick, no direction.
5. 3 tiers protection steel cover and get the ideal centrifugation result.

SPECIFICATIONS

Model	BFH1BY1
Old Model	BCBH-304
Max. Speed	16000rpm
Max. RCF	19040xg
Max. Capacity	6x10ml
Time Range	0~99min59sec
RPM/RCF Convert	Yes
Noise	≤ 65dBA
Speed Accuracy	±20rpm
Acc/Dec	10 Kinds
Size (W x D x Hmm)	370x290x215mm
Temperature accuracy	/
Temperature range	Normal
Net Weight(Kg)	16KG
Power supply	AC 220V/110V 50HZ/60HZ
Certificates	CE,ISO & Calibration report are available
Alt Name	Micro High Speed Centrifuge

ACCESSORIES FOR PURCHASE

No	Name	Max speed (rpm)	Max Volume(ml)	Max RCF (g)
1	Angle rotor	13000	4x8PCR	10410
2	Angle Rotor	16000	12x1.5/2ml	17940
3	Angle Rotor	14000	24x1.5/2ml	15470
4	Angle Rotor	16000	10x5ml	17880
5	Angle Rotor	13000	6x10ml	14190
6	Angle Rotor	16000	40x0.2ml	19040



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com