



LABORATORY LOW SPEED CENTRIFUGE BFH1R1 (BCBLR-405)

LABORATORY LOW SPEED CENTRIFUGE BFH1R1

PRP CENTRIFUGE



1. Frequency conversion motor with simpler construction, more reliable performance, longer life and quietly running.
2. Automatically electric lid lock, super speed, over temperature protection and imbalance protection.
3. Rotor is connected to spindle by specialized taper sleeve, loading simple and quick, no direction, safe and reliable.
4. Microprocessor control, digital display in the speed, time, RCF in operation, speed raising and reducing quick, operate simply.
5. High quality steel body, stainless chamber makes sure the safety of instrument and people.

SPECIFICATIONS

| | |
|--------------------|---|
| Model | BFH1R1 |
| Old Model | BCBLR-405 |
| Max. Speed | 4000rpm |
| Max. RCF | 2820xg |
| Max. Capacity | 4x20ml (PRP tube/kit, brand such as DR.PRP 20ml, Yes/BS prp, 3E prp etc.) |
| Time Range | 0~99min |
| RPM/RCF Convert | Yes |
| Noise (dB) | ≤ 55 |
| Acc/Dec | 10 Kinds |
| Speed Accuracy | ± 20r/min |
| Voltage(V/Hz) | AC 220V/110V 50HZ/60HZ |
| Size (W x D x Hmm) | 540x370x280mm |
| Net Weight(Kg) | 32KG |
| Certificates | CE, ISO & Calibration report are available |
| Alt Name | PRP Centrifuge |



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com