



LABORATORY LOW SPEED CENTRIFUGE BFH1U1 (BCBH-203)

LABORATORY LOW SPEED CENTRIFUGE BFH1U1

TABLETOP HIGH SPEED CENTRIFUGE



1. Brushless DC frequency motor with simpler construction, more reliable performance, longer life and quietly running
2. Automatically electric lid lock, super speed, over temperature protection and imbalance protection. The centrifuge body is made of high quality steel, safe and reliable.
3. Rotor is connected to spindle by specialized taper sleeve, loading simple and quick, no direction, safe and reliable.
4. Microprocessor control.
5. This model with the digital display in the time, speed and RCF.
6. 10 kinds of accelerating and decelerating speed for your choice,
7. 3 layers protection of steel structure, safe and reliable.

SPECIFICATIONS

Model	BFH1U1
Old Model	BCBH-203
Max. Speed	5000rpm
Max. RCF	20600xg
Max. Capacity	6x100ml
Time Range	1~99min
Noise	≤ 60dBA
Speed Accuracy	± 20r/min
Voltage(V/Hz)	AC 220V 50HZ 10A
Size (W x D x Hmm)	513x370x320mm
Net Weight(Kg)	42KG
Package	Carton box
Alt Name	Tabletop High Speed Centrifuge

ACCESSORIES FOR PURCHASE

No	Name	Max speed(rpm)	Max capacity	Max RCF(xg)
1	Angle rotor	16000	12*1.5/2ml	16270
2	Angle rotor	14000	40*0.5ml	19970
3	Angle rotor	15000	24*1.5/2ml	20600
4	Angle rotor	13500	30*1.5/2ml	19340
5	Angle rotor	15000	16x5ml	19350
6	Angle rotor	14000	12x7ml	16370
7	Angle rotor	10000	12x15ml	11840
8	Angle rotor	12000	12x10ml	14510
9	Angle rotor	12000	8x20ml	14510
10	Angle rotor	12000	6x30ml	14000
11	Angle rotor	12000	6x50ml	11140
12	Angle rotor	10000	6x70ml	10810
13	Angle rotor	12000	4x100ml	14850
14	Angle rotor	8000	6x100ml	7280
15	Angle rotor	14000	6x10ml	16460
16	Angle rotor	15000	30x0.5ml	18510
17	Angle rotor	12000	24pieces capillary vessel	15800

18	Angle rotor	11000	48x1.5/2ml	12840
19	Bucket rotor	4000	2x3x48 well	2300



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
 Email: info@biolabscientific.com | Website: www.biolabscientific.com