



MAGNETIC STIRRER BFL1BA1 (BMGS-105)

MAGNETIC STIRRER BFL1BA1

MAGNETIC STIRRER



Speed range of 0-1500 rpm
Maximum stirring quantity of H₂O up to 10L
Ceramic work plate provides excellent corrosion resistance

SPECIFICATIONS

Model	BFL1BA1
Old Model	BMGS-105
Work plate Dimension [W x D]	Ø135mm (5 inch)
Work plate material	Stainless steel cover with ceramic
Motor type	Brushless DC motor
Stirring positions	1
Max. stirring quantity [H ₂ O]	20L
Max. magnetic bar [length]	80mm
Speed range	0-1500rpm
Speed display	Scale
Protection class	IP21
Motor output Power	10W
Power	30W
Voltage, Frequency	100-240V, 50/60Hz
Dimension [WxDxH]	160x280x85mm
Weight	2.8kg
Permissible ambient temperature and humidity	5-40°C, 80%RH
Alt Name	Magnetic Stirrer



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com