



MAGNETIC STIRRER BFL1BA2 (BMGS-108)

MAGNETIC STIRRER BFL1BA2

MAGNETIC STIRRER



Speed range of 0-1500 rpm
Maximum stirring quantity of H₂O up to 10L
Ceramic work plate provides excellent corrosion resistance

SPECIFICATIONS

Model	BFL1BA2
Old Model	BMGS-108
Work plate Dimension [W x D]	184x184mm (7 inch)
Work plate material	Glass ceramic
Motor type	Shaded pole motor
Stirring positions	1
Max. stirring quantity [H ₂ O]	10L
Max. magnetic bar [length]	80mm
Speed range	0-1500rpm
Speed display	Scale
Protection class	IP21
Motor output Power	1.5W
Power	30W
Voltage, Frequency	100-120/200-240V, 50/60Hz
Dimension [WxDxH]	215x360x112mm
Weight	3.8kg
Permissible ambient temperature and humidity	5-40°C, 80%RH
Alt Name	7" Square Magnetic Stirrer



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com