



MAGNETIC STIRRER BFL1BB2 (BMGS-107)

MAGNETIC STIRRER BFL1BB2

MAGNETIC STIRRER



Digital speed control within the range of 0-1500 rpm
Maximum stirring quantity of H₂O at 3L

SPECIFICATIONS

Model	BFL1BB2
Old Model	BMGS-107
Work plate Dimension	Ø135mm (5 inch)
Work plate material	Nylon+GF
Motor type	Brushless DC motor
Stirring positions	1
Max. stirring quantity [H ₂ O]	3L
Max. magnetic bar [length]	50mm
Speed range	0-1500rpm
Speed display	Scale
Protection class	IP42
Motor output Power	3W
Power	10W
Voltage, Frequency	100-240V, 50/60Hz
Dimension [WxDxH]	150x260x80mm
Weight	1.8kg
Permissible ambient temperature and humidity	5-40°C, 80%RH
Alt Name	Magnetic Stirrer



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com