



HOTPLATE MAGNETIC STIRRER BFL1BC1 (BMSS-105)

HOTPLATE MAGNETIC STIRRER BFL1BC1

MAGNETIC HOTPLATE STIRRER



Maintenance-free brushless DC motor
Maximum speed 1100 rpm
Maximum temperature of work plate up to 120°C
Stainless steel work plate, covered with silicone cushion, provides excellent heating uniformity and skid resistance

SPECIFICATIONS

Model	BFL1BC1
Old Model	BMSS-105
Work plate dimension[W x D]	180x450mm
Work Plate material	Stainless steel covered with silicone cushion
Motor type	Brushless DC motor
Stirring positions	10
Max. stirring quantity	0.4Lx10
Max. magnetic bar [length]	40mm
Speed range	0-1100rpm
Speed display	scale
Temperature display	scale
Heating temperature range	Room temp.-120°C
Over temperature protection	140°C
IP Protection class	IP20
Motor rating output	4W
Power	490W
Heating Power	470W
Voltage, Frequency	100-120/200-240V, 50/60Hz
Dimensions[WxDxH]	182x622x65mm
Weight	3.2kg
Permissible ambient temperature and humidity	5-40°C, 80%RH
Alt Name	Magnetic Hotplate Stirrer



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com