



MAGNETIC STIRRER BFL1BD1 (BMSS-104)

MAGNETIC STIRRER BFL1BD1

MAGNETIC STIRRER



Efficient, simple and safe stirring
High performance and space-saving design to maximize the use of benchtop space and saving time by stirring multiple solutions
Maintenance-free brushless DC motor
Speed range 0-1100 rpm
Stainless steel work plate covered with silicone cushion, provides excellent performance of skid resistance

SPECIFICATIONS

Model	BFL1BD1
Old Model	BMSS-104
Work plate Dimension [W x D]	180x450mm
Work plate material	Stainless steel with silicone coating
Motor type	Brushless DC motor
Stirring positions	10
Max. stirring quantity[H2O]	0.4Lx10
Max.magnetic bar[length]	40mm
Speed range	0-1100rpm
Speed display	Scale
Protection class	IP20
Motor output Power	4W
Power	20W
Voltage, Frequency	100-120/200-240V, 50/60Hz
Dimension[WxDxH]	622x205x65mm
Weight	3.2kg
Permissible ambient temperature and humidity	5-40°C, 80%RH
Alt Name	Magnetic Stirrer



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com