



MAGNETIC STIRRER BFL1BE1

MAGNETIC STIRRER BFL1BE1

ULTRA-FLAT COMPACT MAGNETIC STIRRER



Speed range 15 to 1500 rpm
Automatic reverse rotation for excellent mixing
PET top work plate provides excellent chemical resistance
Motorless technology provides ultra silent and vibration-free operation
Perfect magnetic stirrer for use during titration (bottle-top burette) in analytical chemistry

SPECIFICATIONS

Model	BFL1BE1
Work plate Dimension	Ø90mm
Work plate material	PET
Stirring positions	1
Max. stirring quantity[H ₂ O]	0.8L
Max.magnetic bar[length]	40mm
Speed range	15-1500rpm
Protection class	IP21
Output Power	3W
Power	5W
Voltage, Frequency	100-240V, 50/60Hz
Dimension[WxDxH]	175x120x15mm
Weight	0.3kg
Permissible ambient temperature and humidity	5-40°C, 80%RH
Alt Name	Ultra-flat Compact Magnetic Stirrer



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com