



MICRO CENTRIFUGE BFL1K1 (BCMI-401)

MICRO CENTRIFUGE BFL1K1

PALM MICRO CENTRIFUGE



Easy to operate with dual start/stop function - Close lid for a spin speed up to 7000 rpm
Smooth and quiet operation (noise level ≤ 45 dB)
Accelerates and brakes in seconds
Safe, durable and robust rotor with clamp locking design
Compatible with 2 mLx8 and PCR 8x4 rotors
Quick-change rotor without any tool

SPECIFICATIONS

Model	BFL1K1
Old Model	BCMI-401
Max. Speed	7000rpm
Max. RCF	2680xg
Rotor Capacity	0.2/0.5/1.5/2mL x 80.2mL x 32 PCR tubes0.2mL x 4 PCR8-strips
Run Time	Continuous operation
Driving Motor	DC motor
Power Requirements	AC100-240V/50Hz/60Hz20W
Noise Level	≤ 45 dB
Dimension [DxWxH]	160x170x122mm
Weight	0.5kg
Alt Name	Palm Micro Centrifuge

ACCESSORIES FOR PURCHASE

No	Name	Description
1	Rotors	Max. speed: 7000 rpm Max. RCF: 2680 x g Rotor capacity: 0.2 mL x 32 PCR tubes 0.2 mL x 4 PCR 8-strips Rotor material: High strength plastic
2	Rotors	Max. speed: 7000 rpm Max. RCF: 2680 x g Rotor capacity: 2 mL / 1.5 mL x 8 0.2 mL x 8 0.5 mL x 8 Rotor material: High strength plastic
3	Adapters	0.2 mL rotor adapter Used with A8-2P rotor 8 pcs/pk
4	Adapters	0.5 mL rotor adapter Used with A8-2P rotor 8 pcs/pk



1



2



3



4

APPLICATIONS

This centrifuge is ideal for quick spins, microfiltration of samples, cell separations and many other routine laboratory procedures.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com