



## TUBE ROLLER / ROTATOR BFL1T1 (BTUR-201)

# TUBE ROLLER / ROTATOR BFL1T1

## CLASSIC ROTATOR



Speed controlled by Scale  
Adjustable speed range from 0 to 80 rpm

## SPECIFICATIONS

Model	BFL1T1
Old Model	BTUR-201
Motor type	DC motor
Speed range	0-80rpm
Speed display	Scale
Operation mode	Continuous
Voltage, Frequency	100-240V, 50/60Hz
Power	20W
Dimension [DxWxH]	300 x 220 x 310mm
Permissible ambient temperature and humidity	5-40°C, 80%RH
Protection class	IP21
Alt Name	Effective affordable mixing

## ACCESSORIES FOR PURCHASE

No	Name	Description
1	Plastic Clips	Disk/Rotisserie accessory, for 50mL, 15mL and 1.5mL



1

## APPLICATIONS

Laboratory Tube Rotators are ideal for Biochemistry, Molecular Biology and Clinical applications. Provides a gentle but effective mixing of biological samples in 1.5mL to 50mL micro tubes, used in a variety of applications including immune precipitations, prevention of blood coagulation, latex diagnostics etc.



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada  
Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)