



TUBE ROLLER / ROTATOR BFL1W1

TUBE ROLLER / ROTATOR BFL1W1


LCD DIGITAL TUBE ROTATOR



LCD display of speed and time
Circular rotating action, gives gentle but effective mixing
Adjustable speed range from 10 to 70 rpm
Adjustable angle from 0 to 90°
Brushless motor ideal for extended use, and longer life
Wide variety of accessories are available

SPECIFICATIONS

Model	BFL1W1
Motor type	Brushless DC motor
Tilt angle	0-90°
Speed range	10-70rpm
Speed display	LCD
Timer display	LCD
Timer	1-1199min
Operation mode	Continuous / timer
Voltage, Frequency	100-240V, 50/60Hz
Power	40W
Dimension [DxWxH]	280 x 210 x 300mm
Permissible ambient temperature and humidity	5-40°C, 80%RH
Protection class	IP21
Alt Name	LCD Digital Tube Rotator



ACCESSORIES FOR PURCHASE

No	Name	Description
1	Disk Accessories Type-1	for 1.5mLx60 centrifuge tubes holder
2	Disk Accessories Type-2	for 15mLx16 centrifuge tubes holder
3	Disk Accessories Type-3	for 50mLx8 centrifuge tubes holder



1



2



3

APPLICATIONS

Laboratory Tube Rotators are ideal for Biochemistry, Molecular Biology and Clinical applications. Provides a gentle but effective mixing of biological samples in 1.5mL to 50mL micro tubes, used in a variety of applications including immune precipitations, prevention of blood coagulation, latex diagnostics etc.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com