



TYPE I ULTRAPURE WATER PURIFICATION SYSTEM BFN-11001

TYPE I ULTRAPURE WATER PURIFICATION SYSTEM BFN-11001

ULTRAPURE WATER PURIFICATION SYSTEM




This system is a compact solution that combines the production of pure water and ultrapure water in a single unit - eliminating the need for a pretreatment stage upstream of the ultrapure water system.

SPECIFICATIONS

Model	BFN-11001
Production Rate	30L/hr
Water Quality	Type I, Type III
Feed Water	Tap Water
Tank Capacity	30L
Feed Water Requirement :	
Source	Tap water
Conductivity	<2000 $\mu\text{S/cm}$
Hardness	<450 ppm as CaCO_3
Pressure	0.05-0.5 MPa (7-72 psi)
Temperature	5-40°C
Purification Water (Class III) :	
Ionic Rejection	>95%
Bacteria Rejection	>99%
Conductivity Class III	1-20 $\mu\text{S/cm}$ (RO-2: 5 $\mu\text{S/cm}$)
Productivity Rate	30 L/h
Ultrapurification Water (Class I) :	
Resistivity at 25°C	18.2 M $\Omega\cdot\text{cm}$
Conductivity at 25°C	0.055 $\mu\text{S/cm}$
TOC Level	1-5 ppb
Endotoxin (Pyrogens)	N/A
Particulate ($\geq 0.02 \mu\text{m}$)	<1 pc/ml
Bacteria	<0.1 cfu/ml
RNase / DNase	N/A
Manual Dispense Flow Rate	1.5-2.0 L/min
Automatic Dispense Volume	100-60000 ml
Electrical Requirements :	
Electrical Voltage	110 V/220 V $\pm 10\%$
Electrical Frequency	50/60 Hz
Packing Information :	
Net Weight - Main Unit	24 kg
Net Weight - Water Tank (30 L)	7 kg

External Dimensions (W x D x H) Main Unit	315 x 525 x 570 mm
External Dimensions (W x D x H) Water Tank (30 L)	380 x 380 x 595 mm
Shipping Weight Main Unit	37 kg
Shipping Weight Water Tank (30 L)	15 kg
Shipping Dimensions (W x D x H) Main Unit	525 x 610 x 770 mm
Shipping Dimensions (W x D x H) Water Tank (30 L)	520 x 440 x 615 mm
Notes	*If feed water quality is poor (Conductivity > 1000 $\mu\text{S}/\text{cm}$), strengthened pretreatment module and RO-2 type are highly recommended.**When feed water hardness is high (> 450 ppm as CaCO_3), 0.5T water softener is recommended.***Dual-wave UV module needs to be adopted.Also dependent on feed water, recommended feed water TOC < 30 ppb.****The Ultrafiltration module needs to be adopted.Feed water needs to be satisfied as above.
Alt Name	Ultrapure Water Purification System



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com