



HOT AIR STERILIZER BFO1C1 (BSHA-209)

HOT AIR STERILIZER BFO1C1

PORTABLE PRESSURE STEAM STERILIZER

Durable, corrosion resistant and easy to clean stainless steel design in a wide range of capacities. Our hot air sterilizers are microprocessor controlled and perform reliably. A digital display with a broad observation window allows easy operation.



SPECIFICATIONS

| | |
|------------------------------|------------------------------------|
| Model | BFO1C1 |
| Old Model | BSHA-209 |
| Capacity | 18 L |
| Function | Electric heated |
| Sterilization Temperature | 126°C |
| Steam Exhausting Temperature | ≤1°C |
| Working Pressure | 0.14-0.16Mpa |
| Maximum Working Pressure | 0.165 Mpa |
| Chamber Dimension | 280x260 mm |
| Overall Dimension (LxWxH) | 410x410x520 mm |
| Chamber Material | Stainless Steel |
| Weight | 14/16 kg |
| Power | 2 KW |
| Power Supply | 220 V, 50 Hz |
| Alt Name | Portable Pressure Steam Sterilizer |

APPLICATIONS

Medical, Research Institutions, Laboratory, Drinking Water Plant, Food industry, Hospital, Agriculture, Hospital, Scientific Research.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com