



## HOT AIR STERILIZER BFO1C2 (BSHA-210)

# HOT AIR STERILIZER BFO1C2

## PORTABLE PRESSURE STEAM STERILIZER

Durable, corrosion resistant and easy to clean stainless steel design in a wide range of capacities. Our hot air sterilizers are microprocessor controlled and perform reliably. A digital display with a broad observation window allows easy operation.



## SPECIFICATIONS

Model	BFO1C2
Old Model	BSHA-210
Capacity	24 L
Function	Electric heated
Sterilization Temperature	126°C
Steam Exhausting Temperature	≤1°C
Working Pressure	0.14-0.16Mpa
Maximum Working Pressure	0.165 Mpa
Chamber Dimension	280x390mm
Overall Dimension (LxWxH)	410x410x630 mm
Chamber Material	Stainless Steel
Weight	15/17 kg
Power	2 KW
Power Supply	220 V, 50 Hz
Alt Name	Portable Pressure Steam Sterilizer

## APPLICATIONS

Medical, Research Institutions, Laboratory, Drinking Water Plant, Food industry, Hospital, Agriculture, Hospital, Scientific Research.



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)