



ULTRASONIC HOMOGENIZER BFR1C2 (BLIH-402)

ULTRASONIC HOMOGENIZER BFR1C2

PORTABLE ULTRASONIC PROCESSOR

The portable Ultrasonic Processor is used to emulsify, disperse, homogenization, extraction, anti-foaming, cleaning, preparation of nano-materials, dispersion, and speed up chemical reactions.



SPECIFICATIONS

Model	BFR1C2
Old Model	BLIH-402
Operating frequency range	20-25 KHz, frequency automatic tracking
The end of the horn diameter	Φ 3 or 5
Ultrasonic Power	250 W, can be adjusted, fine-tune in steps of 5 W
Parameters sample handling capacity	0.1-100 ml (need to select the appropriate probe)
Generator output power	250 W (max)
Operating voltage	220 V ±5%
Alt Name	Portable Ultrasonic Processor



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com