



## -86°C ULT FREEZER BFRZ-86-202

# -86°C ULT FREEZER BFRZ-86-202

## -86°C ULT CHEST FREEZER



Accurate temperature controller, setting temperature range is -40 to -86°C, accuracy is 0.1 °C.

Suggested ambient temperature: +10 to 32°C.

Stainless Steel inner chamber for friendly use and easy cleaning.

Brand compressor to ensure excellent stability of the cooling system.

Patented refrigeration technology, high-efficiency mixed gas for the cooling system.

Security

Safety door lock to prevent unauthorized access.

Visual and audible alarm system: High/low temperature, sensor failure, controller error.

Power failure protection: turn-on delay of the cooling system after power failure

## SPECIFICATIONS

Model	BFRZ-86-202
Temperature Range	-40 ~ -86°C
Ambient Temperature	10~32°C
Capacity (L/Cu.Ft.)	485 / 17.13
Interior Dimensions (WxDxH) mm	1455x525x674
Exterior Dimensions (WxDxH) mm	1695x840x962
Weight	130 kgs
Noise	55 dB
Cold Rack Max. Loading Capacity (pcs)	2B1U Rack / 24, 2B1D Rack / 6
2" Cold Box Max. Loading Capacity (pcs)	300
0.5~2ml Vials Max. Loading Capacity (pcs)	30,000
Controller	Microprocessor
Display	Digital
Sensor	PTC
Voltage and Frequency	220V±10%, 50/60Hz
Power	670 W
Electric Current	3.16 A
Power Consumption (kWh/24h)	13.5
High / Low Temp. Alarm	Yes
Sensor Error Alarm	Yes
Controller Error Alarm	Yes
Refrigeration Type	Direct Cooling
Defrost	Manual
Compressor (Brand/pcs)	SECOP / 1
Refrigerant / Weight (kgs)	R4001 / 340g
Door (pcs)	Foamed Solid Door / 1
Interior Cabinet Material	Stainless Steel
Exterior Cabinet Material	Coated Steel
Insulation Material	PURF (≥100mm)
Door Lock	Yes
Test Hole (pcs)	1
Fixed / Levelling Castors (pcs)	4
Temperature Datalogger	Optional

Protective Gloves	Optional
Package	Carton box with pallet
Shipping Weight	150 kgs
Shipping Size (WxDxH)	1770x900x1100
Shipping Volume	1.75 CBM
Alt Name	-86°C Ultra-low Freezer



## FEATURES

Accurate temperature controller, setting temperature range is -40 to -86°C, accuracy is 0.1 °C.

Suggested ambient temperature: +10 to 32°C.

Stainless Steel inner chamber for friendly use and easy cleaning.

Brand compressor to ensure excellent stability of the cooling system.

Patented refrigeration technology, high-efficiency mixed gas for the cooling system.

### Security

Safety door lock to prevent unauthorized access.

Visual and audible alarm system: High/low temperature, sensor failure, controller error.

Power failure protection: turn-on delay of the cooling system after power failure.

Controller error protection: the cooling system will remember the normal working cycle and keep working depending on this memory when the controller or two sensors failed.

### Humanized design

Casters with stoppers for easy moving and fixation.

Special-designed service-free drainage system.

## APPLICATIONS

The -86°C ultra-low temperature freezers are designed for the long-term preservation of viruses, germs, vaccines, biological tissues and organs, special food, medicines, reagents, etc.



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)