



SAMPLE TUBE SCANNER BFS-4401

SAMPLE TUBE SCANNER BFS-4401

SAMPLE TUBE SCANNER

Sample tube scanner is specifically designed to operate on top of the sample tubes without applying any weight, avoiding any damage.



It can be finished within 5 seconds for one of 96 holes rack or 384 holes rack. Advanced decoder technology can be identified as stable and secure even if the code is fuzzy.

Optional item for 1-D code rack scanning function.

It can scan the difference of height of the sample tube in the same rack .

Compact, ergonomic design assures efficiency of workflow.

Detection of rack type automatically for setting optimum scanning.

SPECIFICATIONS

Model	BFS-4401
Rack type available	13,24,48,96,384 holes rack
1-D Code available	2/5 Industrial / Interleaved, Code 39, Code 128, Pharmacode, Codabar, EAN 13
2-D Code available	Datamatrix ECC 200, PDF417, QR Code
Camera	10 mega-pixel camera CMOS
Alt Name	Sample tube scanner

APPLICATIONS

Laboratory, Research, Medical



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com