



## PIPETTE PUMP BFS1BH2

## PIPETTE PUMP BFS1BH2

### PIPETTE PUMP

We are engaged in offering a wide range of Pipette Pumps that provide safe and accurate pipetting and dispensing of liquids. It can be rotated to dispense and draw up the liquid. Available in a variety of sizes and shapes.



Economical  
Highly gas tightness  
Environment- friendly materials  
Lightweight and dexterous for operation  
Durable  
Recalibration Mechanism

## SPECIFICATIONS

Model	BFS1BH2
Volume	10ml
Colour	Green + White
Unit	1000
Alt Name	Pipette Pump

## APPLICATIONS

Liquid Handling Task, Institutes, Pharmaceutical, Industrial, Microbiology, Chemistry, Biology, Medical, Laboratory, Research



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)