



TIPS BOX BFS1BO4 (BPIC-704)

TIPS BOX BFS1B04

TIPS BOX

It offers quick and easy access to pipette tips, saving valuable time. Space-saving and environmentally compatible. Commonly used in liquid handling and molecular biology.



Made from medical grade Polypropylene
Slippy surface with excellent design
DNase-free
Rnase-free
Non-pyrogenic
Filter and non-filter for choice
Bag and rack are available to pack

SPECIFICATIONS

Model	BFS1B04
Old Model	BPIC-704
Volume	5000 µl
Unit	20
Alt Name	Tips box

APPLICATIONS

Institutes, Pharmaceutical, Industrial, Microbiology, Plasmid Purification, Chemistry, Biology, Medical, Laboratory, Research, Liquid Handling Task



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com