



## TIPS BOX BFS1BP2 (BPIC-707)

## TIPS BOX BFS1BP2

### STACKABLE TIPS BOX

It offers quick and easy access to pipette tips, saving valuable time. Space-saving and environmentally compatible. Commonly used in liquid handling and molecular biology.



Made from medical grade Polypropylene  
Slippy surface with excellent design  
DNase-free  
Rnase-free  
Non-pyrogenic  
Filter and non-filter for choice  
Bag and rack are available to pack

## SPECIFICATIONS

|           |                    |
|-----------|--------------------|
| Model     | BFS1BP2            |
| Old Model | BPIC-707           |
| Volume    | 200 µl             |
| Unit      | 20                 |
| Alt Name  | Stackable Tips box |

## APPLICATIONS

Institutes, Pharmaceutical, Industrial, Microbiology, Plasmid Purification, Chemistry, Biology, Medical, Laboratory, Research, Liquid Handling Task



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)