




TIPS BOX BFS1BQ6 (BPIC-717)

TIPS BOX BFS1BQ6

TIP BOX WITH ULTRA LOW RETENTION TIP

It offers quick and easy access to pipette tips, saving valuable time. Space-saving and environmentally compatible. Commonly used in liquid handling and molecular biology.



Made from medical grade Polypropylene
Slippy surface with excellent design
DNase-free
Rnase-free
Non-pyrogenic
Filter and non-filter for choice
Bag and rack are available to pack

SPECIFICATIONS

Model	BFS1BQ6
Old Model	BPIC-717
Volume	10000 µl
Unit	10
Alt Name	Tip box with ultra low retention tip

APPLICATIONS

Institutes, Pharmaceutical, Industrial, Microbiology, Plasmid Purification, Chemistry, Biology, Medical, Laboratory, Research, Liquid Handling Task



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com