



## PIPETTE STAND BFS1BU1 (BPIC-302)

## PIPETTE STAND BFS1BU1

### PIPETTE STAND

Pipette stand is a holder manufactured to provide safe and convenient storage for pipettes. Modern design enables this stand to securely hold virtually every pipettes in the marketplace, making it truly 'universal'.



Linear type

Robust design

The freely rotating stand provides safe storage for all micropipettes

The stands are convenient to hold up to 6 pipettes both of single and multi-channel.

## SPECIFICATIONS

Model	BFS1BU1
Old Model	BPIC-302
Description	For 6 pcs PS2 Series Fixed pipettes
Alt Name	Universal Pipette Stand

## APPLICATIONS

Chemistry, Biology, Medical, Laboratory, Research, Liquid Handling Task, Institutes, Pharmaceutical, Industrial, Microbiology, Plasmid Purification



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)