



## TUBE ROLLER / ROTATOR BFS1G1 (BTUR-401)

# TUBE ROLLER / ROTATOR BFS1G1

## DIGITAL ROLLER MIXER

Digital Roller Mixer provides a gentle but highly efficient rocking and rolling motions, ideal for mixing blood samples, viscous substances and liquid-solid suspensions. Used in a variety of applications, prevention of blood coagulation and immune precipitation etc.



Rocking and rolling motions for complete mixing  
Display shows the speed and time  
Speed range of 10-99 rpm  
Seven-roller design with small footprint and easy cleaning  
Digital setting of speed and time

## SPECIFICATIONS

Model	BFS1G1
Old Model	BTUR-401
Speed	Adjustable, 10-99 rpm ( $\pm 10\%$ )
Motion	Simultaneous roll & tilt ( $\pm 5\%$ )
Operating temperature	4-40 °C
Electrical	80-230 V, 50-60 Hz, 50 W
MOQ-10PCS	
Alt Name	Digital Roller Mixer

## FEATURES

Rocking and rolling motions for complete mixing  
Display shows the speed and time  
Speed range of 10-99 rpm  
Seven-roller design with small footprint and easy cleaning  
Digital setting of speed and time



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)