



AIR-JACKETED MULTI GAS INCUBATOR BFS111 (BIMA-113)

AIR-JACKETED MULTI GAS INCUBATOR BFS111

BIO-INCUBATOR

The Biological Indicator Incubator is compact, versatile, reliable and easy to use. Designed for single temperature incubation of the EO Biological Indicators. The thermometer is included with the incubator, providing easy verification of incubation temperature.



Compact, fits in the palm of your hand
LCD display .it s easy to setup and use
Different kinds of aluminum blocks are available
Clear cover ensures temperature uniformity
Portable;Optional 24V adapter

SPECIFICATIONS

| | |
|------------------------------|---------------------------------|
| Model | BFS111 |
| Old Model | BIMA-113 |
| Operating temperature range | From 5°C above ambient to 100°C |
| Time range | 1 min - 99 h 59 min |
| Temp. control accuracy | $\pm 0.5^{\circ}\text{C}$ |
| Display accuracy | 0.1°C |
| Heating time | < 6 min (from 20°C to 57°C) |
| Temp. control accuracy @57°C | $\pm 0.35^{\circ}\text{C}$ |
| Dimensions (WxDxH, mm) | 120 x 152 x 112 |
| MOQ-5PCS | |
| Alt Name | Bio-Incubator |

ACCESSORIES FOR PURCHASE

| No | Name | Specification |
|----|------------|-----------------------------------|
| 1 | tube block | Ø10 mm x 32 mm x 15 (MOQ-5PCS) |
| 2 | tube block | Ø12 mm x 32 mm x 12 (MOQ-5PCS) |
| 3 | tube block | Ø15 mm x 32 mm x 8 (MOQ-5PCS) |
| 4 | tube block | Ø16 mm x 32 mm x 8 (MOQ-5PCS) |
| 5 | tube block | Ø9 mm x 18 (MOQ-5PCS) |



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com