



FLOOR STANDING SHAKING INCUBATOR BFS1N1 (BSTH-109)

FLOOR STANDING SHAKING INCUBATOR BFS1N1

MICROPLATE SHAKER

Shaker for microplates is a type of high speed shaker, which is with brushless DC motor and microprocessor controlled technique, specially for Elisa plates(96/384 wells, tissue culture plates (24/48/96 wells). The shaker is easy to operate with one touch knob and it can load 4 microplates.



Continuous or timed operation
Easy to operate with one touch knob
Variable speed between 200 to 1,500 rpm
DC brushless motor for gentle and powerful motion
Orbital of 3 mm for efficient mixing in small sample vessels

SPECIFICATIONS

Model	BFS1N1
Old Model	BSTH-109
Speed Range	200-1500 RPM
Timing Range	1 min ~ 99 h 59 min
Speed Increment	10 RPM
Mixing Orbit	3 mm
Capacity	4 microplate
Power	220 V / 120 V, 2 A
Dimensions (L x W x H) mm	280 x 270 x 110 mm
Weight	6.5 kg
MOQ-5PCS	
Alt Name	Microplate Shaker



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com