



SINGLE CHANNEL PIPETTE BFS1U2 (BPIP-541)

SINGLE CHANNEL PIPETTE BFS1U2

FIXED VOLUME MICROPIPETTES

New design for smooth activation, excellent ergonomics, easy maintenance and fully autoclavable.



Precision calibration
Recalibration mechanism
Calibration adjustment with simple tools!

SPECIFICATIONS

Model	BFS1U2
Old Model	BPIP-541
Volume	10 μ l
Accuracy $\leq$ %	1.0%
Accuracy $\pm$ μ l	0.1
CV $\leq$ %	0.50%
CV $\pm$ μ l	0.05
MOQ-20PCS	
Alt Name	Fixed Volume Single Channel Micropipette (Ejector 121°C autoclavable)

FEATURES

Precision calibration
Recalibration mechanism



Calibration adjustment with simple tools!



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com