



SINGLE CHANNEL PIPETTE BFS1U4 (BPIP-543)

SINGLE CHANNEL PIPETTE BFS1U4

FIXED VOLUME MICROPIPETTES

New design for smooth activation, excellent ergonomics, easy maintenance and fully autoclavable.



Precision calibration
Recalibration mechanism
Calibration adjustment with simple tools!

SPECIFICATIONS

Model	BFS1U4
Old Model	BIPI-543
Volume	25 µl
Accuracy ≤ %	1.0%
Accuracy ± µl	0.25
CV ≤ %	0.30%
CV ± µl	0.075
MOQ-20PCS	
Alt Name	Fixed Volume Single Channel Micropipette (Ejector 121°C autoclavable)

FEATURES

Precision calibration
Recalibration mechanism



Calibration adjustment with simple tools!



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com