



## ANALYTICAL BALANCE BFT-1102

# ANALYTICAL BALANCE BFT-1102

## INTERNAL CALIBRATION ANALYTICAL BALANCE



Internal calibration  
LCD (White back light with black font)  
One key / Time trigger internal calibration  
Electromagnetic force sensor  
Glass draft shield with 3 sliding doors for easy access to the items being weighed  
Tare function / Counting / Percentage  
Unit conversion (g/mg/ct/oz...)  
Sensitivity / Speed set  
Overload alarm / Fault alarm / Level indicator

## SPECIFICATIONS

|                   |                                         |
|-------------------|-----------------------------------------|
| Model             | BFT-1102                                |
| Capacity          | 410g                                    |
| Readability       | 0.1mg                                   |
| Min weighing      | 0.4mg                                   |
| Stable Time       | ≤4s                                     |
| Display           | LCD(white back light with black font)   |
| Pan size          | Φ80mm                                   |
| Operation temp    | 5-35°C                                  |
| Repeat ability    | 0.1mg if weight ≤ 220g                  |
| Linearity         | ± 0.5mg                                 |
| Draft shield size | 162mmx171mmx225mm                       |
| Cal.weight        | Internal calibration                    |
| Interface         | RS232(Standard)                         |
| Dimension         | 345x215x345mm                           |
| N. weight         | 7 kg                                    |
| Alt Name          | Internal Calibration Analytical Balance |

## ACCESSORIES FOR PURCHASE

| No | Name       | Description                                             | Information                                                                                                                                      |
|----|------------|---------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------|
| 1  | Weight set | Recommend corresponding accuracy calibration (1mg-20kg) | M1 levelSuitable for 0.1g-0.01g<br>F2 levelSuitable for 0.01g-0.001g<br>F1 levelSuitable for 0.001g-0.0001g<br>E2 levelSuitable for over 0.0001g |
| 2  | Weight set | Lock type series weights (5kg-20kg)                     | M1 levelSuitable for 0.1g-0.01g<br>F2 levelSuitable for 0.01g-0.001g<br>F1 levelSuitable for 0.001g-0.0001g                                      |



1



1




1



1



2



2



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)