



## ANALYTICAL BALANCE BFT-1507

# ANALYTICAL BALANCE BFT-1507

## ANALYTICAL BALANCE



Electromagnetic force sensor  
Aluminum shell / Large space wind proof cover  
LCD (White back light with black font)  
Tare function / Counting / Unit conversion (g/mg/ct/oz...)  
Sensitivity / Speed set / Percentage  
Overload alarm / Fault alarm / Level indicator  
Option: Under weighing / Printer

## SPECIFICATIONS

|                   |                                            |
|-------------------|--------------------------------------------|
| Model             | BFT-1507                                   |
| Capacity          | 300g                                       |
| Readability       | 0.0001g                                    |
| Min Weighing      | 0.0004g                                    |
| Repeatability     | ±0.0002g                                   |
| Linearity         | ±0.0003g                                   |
| Stable Time       | ≤3s                                        |
| Pan Size          | Φ80mm                                      |
| Appearance Size   | 365*210*355mm                              |
| Draft Shield Size | 175*165*230mm (Inner diameter)             |
| Communication     | RS232 (Standard) USB (Option)              |
| Temperature       | 15-35°C                                    |
| Calibration       | External calibration                       |
| Cal. Weight       | 200g                                       |
| Option            | Printer/Under weighing                     |
| Packing Size      | Inside:440*295*430mm Outside:490*335*485mm |
| Gross Weight      | 10kg                                       |
| Alt Name          | Analytical Balance                         |

## ACCESSORIES FOR PURCHASE

| No | Name       | Description                                             | Information                                                                                                                                      |
|----|------------|---------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------|
| 1  | Weight set | Recommend corresponding accuracy calibration (1mg-20kg) | M1 levelSuitable for 0.1g-0.01g<br>F2 levelSuitable for 0.01g-0.001g<br>F1 levelSuitable for 0.001g-0.0001g<br>E2 levelSuitable for over 0.0001g |
| 2  | Weight set | Lock type series weights (5kg-20kg)                     | M1 levelSuitable for 0.1g-0.01g<br>F2 levelSuitable for 0.01g-0.001g<br>F1 levelSuitable for 0.001g-0.0001g                                      |



1



1




1



1



2



2



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)